

Wall mounted AED support Automated External Defibrillator

AIVIA «S»

Presentation

The AIVIA «S» is designed for an indoor positioning.

Combining strength and lightness, the AIVIA «S» is composed of 2 distinct parts:


- the polycarbonate cover,
- the ABS chassis.

These materials are well known in the industry, giving excellent shock resistance properties.

Access to the AED is protected by coated seals.

Markings

On the front cover are shown the survival sequence and three clear pictograms with instruction for use.



AIVIA «S»

Dimensions :

H = 400 mm, L = 341 mm, P = 198 mm,

Weight = 1,5 kg



AED holding strap



Protection seals

Wall mount holes 8,5 mm



PYRES.COM

Mas des Tilleuls
66680 Canohès - FRANCE
Tél : +33 (0)4 68 68 39 68
Fax : +33 (0)4 68 68 39 69
aivia@pyres.com
www.pyres.com

Product distributed by :

AIVIA «S»